












Title: Qualitative Organic Chemistry.

Aim: To identify the elements present in an organic compound and determine its functional group.






Experiment	Observation	Inference
1. <u>Soda lime test:</u> Heat substance and soda lime (1:2) in a test tube.	Smell of NH_3 : (a) turns red litmus blue (b) turns mercurous nitrate paper black.	Nitrogen is present.
2. <u>Lassaigne's Test:</u> Fuse a little amount of substance with a pea-size metallic solution in a fusion tube. Heat gently at first, then strongly. Extract the melt with distilled water and filter.		




Experiment	Observation	Inference
a) To 5 ml of filtrate, add 2-3 drops of NaOH, 1-2 drops of FeSO_4 and boil. Add 2-3 drops of FeCl_3 and acidify with HCl.	Blue/green precipitate or colouration (prussian blue)	Nitrogen present.
b) To a part of fit filtrate, add fresh solution of sodium nitroprusside.	Violet precipitate	Sulphur present
To another part, acidify with acetic acid and add lead acetate solution.	Black precipitate	Sulphur present
c) To a 2 ml solution of filtrate, add dil. HNO_3 and boil for a few minutes. Cool and add AgNO_3 solution.	(i) Curdy white precipitate (ii) Curdy yellow precipitate	(i) Chlorine present (ii) Bromine/Iodine present.
d) To a 2 ml of filtrate add H_2SO_4 , then few drops of fresh chlorine water and 2 ml of CHCl_3 or CCl_4 or CS_2 and shake well.	(i) Organic layer turns yellow violet. (ii) Organic layer turns yellow	Iodine present Bromine present
If the organic layer is violet, add more chlorine water and shake well until colour is discharged	Organic layer turns yellow	Iodine and Bromine both present



Experiment	Observation	Inference
3. <u>Test for -COOH group:</u>		
a) Test a little of the sample with blue litmus	Blue litmus turns red	-COOH group present
b) Dissolve or suspend the sample in water, add solid or aqueous solution of NaHCO_3 , small quantity at a time.	Effervescence with evolution of CO_2	
4. <u>Test for Phenolic -OH group:</u>	(i) Reddish colouration or precipitate.	(i) Simple carboxylic acid
Dissolve the substance in water or alcohol and add a drop of FeCl_3 solution	(ii) Green, blue or violet colouration	(ii) Phenolic -OH group present
5. <u>Test for Ester group:</u>		
Add two beads of NaOH to a little of the sample in water minimum volume of water and boil for a few minutes. Cool and then acidify with dilute HCl .	White solid	-COOR group present.



	Experiment	Observation	Inference
6.	<p><u>Test for Carbonyl group:</u></p> <p>2,4-DNP test - To a small portion of alcoholic solution of the substance or to a small amount of the sample, add minimum quantity of glacial acetic acid, warm till clear solution appears. Then add 2 ml of standard solution of 2,4-dinitrophenyl hydrazine. Shake gently and heat to boiling for few minutes. Allow to cool. If no precipitate appears, add 2 drops of conc. H_2SO_4. Warm for about five minutes again and scratch inside of the test tube with a glass rod. Cool and allow to stand.</p>	<p>Orange or red precipitate.</p>	<p>Carbonyl group present ($>C=O$ or $-CHO$)</p>



Experiment	Observation	Inference
7. <u>Test for primary amine:</u> Diazotisation - Add a pinch of NaNO_2 to a well cooled HCl solution of the sample. Then pour a little of diazotised solution into excess of cold alkaline solution of β -naphthol.	Immediately brilliant scarlet or red precipitate appears.	$-\text{NH}_2$ group present
8. <u>Test for amide group:</u> (i) Heat 0.5 g of the sample with 2 ml 50% NaOH solution	(i) Smell of NH_3 - which turns $\text{Hg}(\text{NO}_3)_2$ paper black black	$-\text{CONH}_2$ group present
(ii) Treat a little of the aqueous solution of the sample with a few drops of cold solution of HNO_2 (HCl and NaNO_2).	(ii) N_2 gas evolves	

